

Anti-Human/Mouse VPS26B Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA06128

3th Edition

Description: VPS26B (VPS26, Retromer Complex Component B) is a Protein Coding gene. Among its related pathways are Endocytosis. GO annotations related to this gene include protein transporter activity. An important paralog of this gene is VPS26A.

Antigen: Synthesized peptide derived from the Internal region of human VPS26B

Form:

How to use: 1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution: PBS (pH7.4) containing 1% BSA

Application: This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: