

Anti-Human VRK3 Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA06129

3th Edition

Description: Inactive kinase that suppresses ERK activity by promoting phosphatase activity of DUSP3 which specifically dephosphorylates and inactivates ERK in the nucleus.

Antigen: Synthesized peptide derived from the Internal region of human VRK3.

Form:

How to use: 1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution: PBS (pH7.4) containing 1% BSA

Application: This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: