

**Anti-Human/Mouse/Rat VTI1B Polyclonal Antibody**

**Polyclonal Antibody**

**Cat.NO.: PA06130**

---

3th Edition

**Description:**V-SNARE that mediates vesicle transport pathways through interactions with t-SNAREs on the target membrane. These interactions are proposed to mediate aspects of the specificity of vesicle trafficking and to promote fusion of the lipid bilayers. May be concerned with increased secretion of cytokines associated with cellular senescence.

**Antigen:**Synthesized peptide derived from the Internal region of human v-SNARE Vti1p.

**Form:**

**How to use:**1.0 ml distilled water will be added to the product

**Stability:** Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

**Dilution:**PBS (pH7.4) containing 1% BSA

**Application:**This antibody can be used for western blotting in concentration of 1?5?g/ml.

**Specificity:**Expressed in all tissues examined.