

Anti-Human/Mouse ZAP70 Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA06155

3th Edition

Description: Plays a role in T-cell development and lymphocyte activation. Essential for TCR-mediated IL-2 production. Isoform 1 induces TCR-mediated signal transduction, isoform 2 does not.

Antigen: Synthesized peptide derived from human ZAP-70 around the non-phosphorylation site of Y292.

Form:

How to use: 1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution: PBS (pH7.4) containing 1% BSA

Application: This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: Expressed in T- and natural killer cells.