

Anti-Human ZAR1 Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA06156

3th Edition

Description: This maternal effect gene is oocyte-specific and encodes a protein that is thought to function in the initiation of embryogenesis. A similar protein in mouse is required for female fertility. ZAR1 (Zygote Arrest 1) is a Protein Coding gene. Among its related pathways are Preimplantation Embryo. An important paralog of this gene is ZAR1L.

Antigen: Synthesized peptide derived from the Internal region of human ZAR1.

Form:

How to use: 1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution: PBS (pH7.4) containing 1% BSA

Application: This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: Ovary and testis.