

**Anti-Human/Monkey ZC3H11A Polyclonal Antibody**

**Polyclonal Antibody**

**Cat.NO.: PA06159**

---

3th Edition

**Description:**ZC3H11A (Zinc Finger CCCH-Type Containing 11A) is a Protein Coding gene. Among its related pathways are Gene Expression and Transport of Mature Transcript to Cytoplasm. GO annotations related to this gene include poly(A) RNA binding. An important paralog of this gene is ZC3H11B.

**Antigen:**Synthesized peptide derived from the Internal region of human ZC3H11A

**Form:**

**How to use:**1.0 ml distilled water will be added to the product

**Stability:** Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

**Dilution:**PBS (pH7.4) containing 1% BSA

**Application:**This antibody can be used for western blotting in concentration of 1?5?g/ml.

**Specificity:**