

Anti-Human ZNF397 Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA06172

3th Edition

Description: Isoform 3 acts as a DNA-dependent transcriptional repressor.

Antigen: Synthesized peptide derived from the N-terminal region of human ZNF397

Form:

How to use: 1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution: PBS (pH7.4) containing 1% BSA

Application: This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: Expressed strongly in testis, moderately in skeletal muscle, pancreas and prostate, and weakly in heart, placenta, liver, kidney, spleen, thymus and small intestine.