

**Anti-Human/Mouse ZNF420 Polyclonal Antibody**

**Polyclonal Antibody**

**Cat.NO.: PA06173**

---

3th Edition

**Description:**ZNF420 contains 19 C2H2-type zinc fingers and 1 KRAB domain. It may be involved in transcriptional regulation. It is a negative regulator of p53-mediated apoptosis.

**Antigen:**Synthesized peptide derived from the Internal region of human ZNF420

**Form:**

**How to use:**1.0 ml distilled water will be added to the product

**Stability:** Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

**Dilution:**PBS (pH7.4) containing 1% BSA

**Application:**This antibody can be used for western blotting in concentration of 1?5?g/ml.

**Specificity:**