

Anti-Human ZNF436 Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA06174

3th Edition

Description: May be a transcriptional repressor.

Antigen: Synthesized peptide derived from the Internal region of human ZNF436

Form:

How to use: 1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution: PBS (pH7.4) containing 1% BSA

Application: This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: Expressed in fetal brain, heart, liver, spleen, bladder, lung, skin, skeletal muscle, stomach and pancreas.