

Anti-Human ZNF460 Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA06175

3th Edition

Description: May be involved in transcriptional regulation.

Antigen: Synthesized peptide derived from the C-terminal region of human ZNF460

Form:

How to use: 1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution: PBS (pH7.4) containing 1% BSA

Application: This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: Ubiquitously expressed at low levels. Highest levels are found in pancreas and liver.