

Anti-Human/Mouse ZNHIT1 Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA06185

3th Edition

Description:ZNHIT1, Zinc finger HIT domain containing protein 1, appears to play a role in p53-mediated apoptosis induction.

Antigen:Synthesized peptide derived from the Internal region of human ZNHIT1

Form:

How to use:1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution:PBS (pH7.4) containing 1% BSA

Application:This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity:Expressed abundantly in liver, but weakly in skeletal muscle, ovary and small intestine.