

Anti-Human/Mouse/Rat ZnT-9 Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA06187

3th Edition

Description: SLC30A9 (Solute Carrier Family 30 Member 9) is a Protein Coding gene. Among its related pathways are Metal ion SLC transporters. GO annotations related to this gene include transcription factor activity, sequence-specific DNA binding and ligand-dependent nuclear receptor transcription coactivator activity.

Antigen: Synthesized peptide derived from the Internal region of human ZnT-9

Form:

How to use: 1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution: PBS (pH7.4) containing 1% BSA

Application: This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: Ubiquitously expressed in fetal and adult tissues and cancer cell lines.