

Anti-Human FAIM2 Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA06230

3th Edition

Description:Protects cells uniquely from Fas-induced apoptosis.

Antigen:Synthesized peptide derived from the N-terminal region of human FAIM2.

Form:

How to use:1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution:PBS (pH7.4) containing 1% BSA

Application:This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity:Widely expressed. Expressed at high levels in the brain especially in the hippocampus.