

**Anti-Human/Mouse/Rat HGFAC Polyclonal Antibody**

**Polyclonal Antibody**

**Cat.NO.: PA06232**

---

3th Edition

**Description:** This gene encodes a member of the peptidase S1 protein family. The encoded protein is first synthesized as an inactive single-chain precursor before being activated to a heterodimeric form by endoproteolytic processing. It acts as serine protease that converts hepatocyte growth factor to the active form. Alternative splicing results in multiple transcript variants.

**Antigen:** Synthesized peptide derived from the C-terminal region of human HGFA.

**Form:**

**How to use:** 1.0 ml distilled water will be added to the product

**Stability:** Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

**Dilution:** PBS (pH7.4) containing 1% BSA

**Application:** This antibody can be used for western blotting in concentration of 1?5?g/ml.

**Specificity:** Liver.