

**Anti-Human/Mouse/Rat CX3CR1 Polyclonal Antibody**

**Polyclonal Antibody**

**Cat.NO.: PA06267**

---

3th Edition

**Description:**Fractalkine is a transmembrane protein and chemokine involved in the adhesion and migration of leukocytes. The protein encoded by this gene is a receptor for fractalkine. The encoded protein also is a coreceptor for HIV-1, and some variations in this gene lead to increased susceptibility to HIV-1 infection and rapid progression to AIDS. Four transcript variants encoding two different isoforms have been found for this gene.

**Antigen:**Synthesized peptide derived from the Internal region of human Fractalkine Receptor

**Form:**

**How to use:**1.0 ml distilled water will be added to the product

**Stability:** Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

**Dilution:**PBS (pH7.4) containing 1% BSA

**Application:**This antibody can be used for western blotting in concentration of 1?5?g/ml.

**Specificity:**Expressed in lymphoid and neural tissues.