

Anti-Human/Mouse/Rat RRP7A Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA06268

3th Edition

Description:RRP7A (Ribosomal RNA Processing 7 Homolog A) is a Protein Coding gene. Among its related pathways are rRNA processing in the nucleus and cytosol and Gene Expression. GO annotations related to this gene include poly(A) RNA binding and nucleotide binding.

Antigen:Synthesized peptide derived from the C-terminal region of human RRP7A.

Form:

How to use:1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution:PBS (pH7.4) containing 1% BSA

Application:This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: