

Anti-Human/Mouse/Rat HSP40-4 Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA06313

3th Edition

Description:Co-chaperone of Hsc70. Seems to play a role in protein import into mitochondria.

Antigen:Synthesized peptide derived from the Internal region of human HSP40-4.

Form:

How to use:1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution:PBS (pH7.4) containing 1% BSA

Application:This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity:Ubiquitous. Isoform 2 is highly expressed in testis and lung, but detected at low levels in thymus, prostate, colon and liver.