

Anti-Human Galectin 3 Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA06318

3th Edition

Description:Galactose-specific lectin which binds IgE. May mediate with the alpha-3, beta-1 integrin the stimulation by CSPG4 of endothelial cells migration. Together with DMBT1, required for terminal differentiation of columnar epithelial cells during early embryogenesis.

Antigen:Synthesized peptide derived from the C-terminal region of human CBP 35.

Form:

How to use:1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution:PBS (pH7.4) containing 1% BSA

Application:This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity:A major expression is found in the colonic epithelium. It is also abundant in the activated macrophages.