

Anti-Human/Rat MDK Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA06326

3th Edition

Description:Has heparin binding activity, and growth promoting activity. Involved in neointima formation after arterial injury, possibly by mediating leukocyte recruitment. Also involved in early fetal adrenal gland development.

Antigen:Synthesized peptide derived from the Internal region of human MK.

Form:

How to use:1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution:PBS (pH7.4) containing 1% BSA

Application:This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity:Expressed in various tumor cell lines. In insulinoma tissue predominantly expressed in precancerous lesions.