

Anti-Human/Mouse/Rat CDC34 Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA06338

3th Edition

Description:The protein encoded by this gene is a member of the ubiquitin-conjugating enzyme family. Ubiquitin-conjugating enzyme catalyzes the covalent attachment of ubiquitin to other proteins. This protein is a part of the large multiprotein complex, which is required for ubiquitin-mediated degradation of cell cycle G1 regulators, and for the initiation of DNA replication.

Antigen:Synthesized peptide derived from the Internal region of human Cdc34.

Form:

How to use:1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution:PBS (pH7.4) containing 1% BSA

Application:This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity:Expressed in testes during spermatogenesis to regulate repression of cAMP-induced transcription.