

**Anti-Human/Rat CRP Polyclonal Antibody**

**Polyclonal Antibody**

**Cat.NO.: PA06345**

---

3th Edition

**Description:** Displays several functions associated with host defense: it promotes agglutination, bacterial capsular swelling, phagocytosis and complement fixation through its calcium-dependent binding to phosphorylcholine. Can interact with DNA and histones and may scavenge nuclear material released from damaged circulating cells

**Antigen:** Synthesized peptide derived from the Internal region of human CRP.

**Form:**

**How to use:** 1.0 ml distilled water will be added to the product

**Stability:** Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

**Dilution:** PBS (pH7.4) containing 1% BSA

**Application:** This antibody can be used for western blotting in concentration of 1?5?g/ml.

**Specificity:** Found in plasma.