

**Anti-Human/Mouse/Rat IL1B Polyclonal Antibody**

**Polyclonal Antibody**

**Cat.NO.: PA06350**

---

3th Edition

**Description:**Produced by activated macrophages, IL-1 stimulates thymocyte proliferation by inducing IL-2 release, B-cell maturation and proliferation, and fibroblast growth factor activity. IL-1 proteins are involved in the inflammatory response, being identified as endogenous pyrogens, and are reported to stimulate the release of prostaglandin and collagenase from synovial cells.

**Antigen:**Synthesized peptide derived from the Internal region of human IL1B.

**Form:**

**How to use:**1.0 ml distilled water will be added to the product

**Stability:** Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

**Dilution:**PBS (pH7.4) containing 1% BSA

**Application:**This antibody can be used for western blotting in concentration of 1?5?g/ml.

**Specificity:**Expressed in activated monocytes/macrophages (at protein level).