

Anti-Human/Mouse/Rat FOXP1 Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA06361

3th Edition

Description: Transcriptional repressor. It plays an important role in the specification and differentiation of lung epithelium. Can act with CTBP1 to synergistically repress transcription but CTPBP1 is not essential. Essential transcriptional regulator of B cell development.

Antigen: Synthesized peptide derived from the C-terminal region of human FOXP1.

Form:

How to use: 1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution: PBS (pH7.4) containing 1% BSA

Application: This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: Isoform 8 is specifically expressed in embryonic stem cells.