

**Anti-Human SERPINA4 Polyclonal Antibody**

**Polyclonal Antibody**

**Cat.NO.: PA06363**

---

3th Edition

**Description:**Inhibits human amidolytic and kininogenase activities of tissue kallikrein. Inhibition is achieved by formation of an equimolar, heat- and SDS-stable complex between the inhibitor and the enzyme, and generation of a small C-terminal fragment of the inhibitor due to cleavage at the reactive site by tissue kallikrein.

**Antigen:**Synthesized peptide derived from the Internal region of human Kallistatin.

**Form:**

**How to use:**1.0 ml distilled water will be added to the product

**Stability:** Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

**Dilution:**PBS (pH7.4) containing 1% BSA

**Application:**This antibody can be used for western blotting in concentration of 1?5?g/ml.

**Specificity:**Expressed by the liver and secreted in plasma.