

**Anti-Human/Mouse/Rat LASP1 Polyclonal Antibody**

**Polyclonal Antibody**

**Cat.NO.: PA06381**

---

3th Edition

**Description:** Plays an important role in the regulation of dynamic actin-based, cytoskeletal activities. Agonist-dependent changes in LASP1 phosphorylation may also serve to regulate actin-associated ion transport activities, not only in the parietal cell but also in certain other F-actin-rich secretory epithelial cell types.

**Antigen:** Synthesized peptide derived from the Internal region of human LASP-1.

**Form:**

**How to use:** 1.0 ml distilled water will be added to the product

**Stability:** Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

**Dilution:** PBS (pH7.4) containing 1% BSA

**Application:** This antibody can be used for western blotting in concentration of 1?5?g/ml.

**Specificity:**