

Anti-Human PSG1 Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA06398

3th Edition

Description: The soluble form from urine retains its specific complement binding activity, but exhibits greatly reduced ability to inhibit MAC assembly on cell membranes.

Antigen: Synthesized peptide derived from the Internal region of human PSG1.

Form:

How to use: 1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution: PBS (pH7.4) containing 1% BSA

Application: This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: