



**Anti-Human/Rat CD107a Polyclonal Antibody**

**Polyclonal Antibody**

**Cat.NO.: PA06404**

---

3th Edition

**Description:**The protein encoded by this gene is a member of a family of membrane glycoproteins. This glycoprotein provides selectins with carbohydrate ligands. It may also play a role in tumor cell metastasis.

**Antigen:**Synthesized peptide derived from the Internal region of human CD107a.

**Form:**

**How to use:**1.0 ml distilled water will be added to the product

**Stability:** Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

**Dilution:**PBS (pH7.4) containing 1% BSA

**Application:**This antibody can be used for western blotting in concentration of 1?5?g/ml.

**Specificity:**