



Anti-Human CK-8 Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA06462

3th Edition

Description: Together with KRT19, helps to link the contractile apparatus to dystrophin at the costameres of striated muscle.

Antigen: Synthesized peptide derived from human Cytokeratin 8 around the non-acetylation site of K483.

Form:

How to use: 1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution: PBS (pH7.4) containing 1% BSA

Application: This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: Observed in muscle fibers accumulating in the costameres of myoplasm at the sarcolemma membrane in structures that contain dystrophin and spectrin. Expressed in gingival mucosa and hard palate of the oral cavity.