

Anti-Human/Mouse/Rat CD53 Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA06517

3th Edition

Description: May be involved in growth regulation in hematopoietic cells.

Antigen: Synthesized peptide derived from the Internal region of human CD53.

Form:

How to use: 1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution: PBS (pH7.4) containing 1% BSA

Application: This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: B-cells, monocytes, macrophages, neutrophils, single (CD4 or CD8) positive thymocytes and peripheral T-cells.