

**Anti-Human CXCR5 Polyclonal Antibody**

**Polyclonal Antibody**

**Cat.NO.: PA06536**

---

3th Edition

**Description:** Cytokine receptor that binds to B-lymphocyte chemoattractant (BLC). Involved in B-cell migration into B-cell follicles of spleen and Peyer patches but not into those of mesenteric or peripheral lymph nodes. May have a regulatory function in Burkitt lymphoma (BL) lymphomagenesis and/or B-cell differentiation.

**Antigen:** Synthesized peptide derived from the N-terminal region of human Blr1.

**Form:**

**How to use:** 1.0 ml distilled water will be added to the product

**Stability:** Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

**Dilution:** PBS (pH7.4) containing 1% BSA

**Application:** This antibody can be used for western blotting in concentration of 1?5?g/ml.

**Specificity:** Expression in mature B-cells and Burkitt lymphoma cells.