

Anti-Human CD300LD Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA06554

3th Edition

Description:CD300LD (CD300 Molecule Like Family Member D) is a Protein Coding gene. Among its related pathways are Immunoregulatory interactions between a Lymphoid and a non-Lymphoid cell and Innate Immune System. An important paralog of this gene is CD300LF.

Antigen:Synthesized peptide derived from the Internal region of human CD300d.

Form:

How to use:1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution:PBS (pH7.4) containing 1% BSA

Application:This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity:Expression seems restricted to cells of myeloid lineage.