

Anti-Human/Mouse/Rat FUNDC1 Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA06575

3th Edition

Description: This gene encodes a protein with a FUN14 superfamily domain. The function of the encoded protein is not known. FUNDC1 (FUN14 Domain Containing 1) is a Protein Coding gene. Diseases associated with FUNDC1 include Hypoxia. Among its related pathways are Mitophagy - animal. An important paralog of this gene is FUNDC2.

Antigen: Synthesized peptide derived from the N-terminal region of human FUNDC1.

Form:

How to use: 1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution: PBS (pH7.4) containing 1% BSA

Application: This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: Widely expressed.