

**Anti-Human CD66e Polyclonal Antibody**

**Polyclonal Antibody**

**Cat.NO.: PA06578**

---

3th Edition

**Description:** Carcinoembryonic antigen-related cell adhesion molecule 5 (CEACAM5) also known as CD66e (Cluster of Differentiation 66e), is a human gene. Members of the CEACAM subfamily, including CEACAM5, belong to the CEA gene family. CD66E is a well known tumor marker and a heavily glycosylated GPI-linked cell surface molecule.

**Antigen:** Recombinant protein of human CEACAM5

**Form:**

**How to use:** 1.0 ml distilled water will be added to the product

**Stability:** Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

**Dilution:** PBS (pH7.4) containing 1% BSA

**Application:** This antibody can be used for western blotting in concentration of 1?5?g/ml.

**Specificity:** Found in adenocarcinomas of endodermally derived digestive system epithelium and fetal colon.