

**Anti-Human/Mouse ADAP1 Polyclonal Antibody**

**Polyclonal Antibody**

**Cat.NO.: PA06605**

---

3th Edition

**Description:**Arf-GAP with dual PH domain-containing protein 1 is a protein that in humans is encoded by the ADAP1 gene.GTPase-activating protein for the ADP ribosylation factor family. Binds phosphatidylinositol 3,4,5-trisphosphate (PtdInsP3) and inositol 1,3,4,5-tetrakisphosphate (InsP4).

**Antigen:**Recombinant protein of human ADAP1

**Form:**

**How to use:**1.0 ml distilled water will be added to the product

**Stability:** Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

**Dilution:**PBS (pH7.4) containing 1% BSA

**Application:**This antibody can be used for western blotting in concentration of 1?5?g/ml.

**Specificity:**Expressed at highest levels in brain and at lower levels in peripheral blood leukocytes.