

Anti-Human/Mouse COX19 Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA06642

3th Edition

Description:COX19 encodes a cytochrome c oxidase (COX)-assembly protein. The *S. cerevisiae* Cox19 protein may play a role in metal transport to the mitochondrial intermembrane space and assembly of complex IV of the mitochondrial respiratory chain.

Antigen:Recombinant protein of human COX19

Form:

How to use:1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution:PBS (pH7.4) containing 1% BSA

Application:This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity:Ubiquitously expressed. Highly expressed in skeletal muscle.