

Anti-Human/Mouse/Rat ANGPT2 Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA06692

3th Edition

Description:The protein encoded by this gene is an antagonist of angiotensin II (ANGPT1) and endothelial TEK tyrosine kinase (TIE-2, TEK). The encoded protein disrupts the vascular remodeling ability of ANGPT1 and may induce endothelial cell apoptosis. Three transcript variants encoding three different isoforms have been found for this gene.

Antigen:Recombinant protein of human ANGPT2

Form:

How to use:1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution:PBS (pH7.4) containing 1% BSA

Application:This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: