

Anti-Human/Mouse TNFRSF12A Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA06821

3th Edition

Description: Tumor necrosis factor receptor superfamily member 12A also known as the Tweak-receptor (TweakR) is a protein that in humans is encoded by the TNFRSF12A gene. Receptor for TNFSF12/TWEAK. Weak inducer of apoptosis in some cell types. Promotes angiogenesis and the proliferation of endothelial cells. May modulate cellular adhesion to matrix proteins. Highly expressed in heart, placenta and kidney. Intermediate expression in lung, skeletal muscle and pancreas.

Antigen: Recombinant protein of human TNFRSF12A

Form:

How to use: 1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution: PBS (pH7.4) containing 1% BSA

Application: This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: Highly expressed in heart, placenta and kidney. Intermediate expression in lung, skeletal muscle and pancreas.