

**Anti-Human/Mouse GAS2 Polyclonal Antibody**

**Polyclonal Antibody**

**Cat.NO.: PA06876**

---

3th Edition

**Description:**The protein encoded by this gene is a caspase-3 substrate that plays a role in regulating microfilament and cell shape changes during apoptosis. It can also modulate cell susceptibility to p53-dependent apoptosis by inhibiting calpain activity. Multiple alternatively spliced variants, encoding the same protein, have been identified.

**Antigen:**Recombinant protein of human GAS2

**Form:**

**How to use:**1.0 ml distilled water will be added to the product

**Stability:** Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

**Dilution:**PBS (pH7.4) containing 1% BSA

**Application:**This antibody can be used for western blotting in concentration of 1?5?g/ml.

**Specificity:**Ubiquitously expressed with highest levels in liver, lung, and kidney. Not found in spleen.