

Anti-Human/Mouse HOXA6 Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA06903

3th Edition

Description:In vertebrates, the genes encoding the class of transcription factors called homeobox genes are found in clusters named A, B, C, and D on four separate chromosomes. Expression of these proteins is spatially and temporally regulated during embryonic development. This gene is part of the A cluster on chromosome 7 and encodes a DNA-binding transcription factor which Human, Mouse regulate gene expression, morphogenesis, and differentiation.

Antigen:Recombinant protein of human HOXA6

Form:

How to use:1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution:PBS (pH7.4) containing 1% BSA

Application:This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: