

Anti-Human/Mouse HSF1 Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA06908

3th Edition

Description:The product of this gene is a heat-shock transcription factor. Transcription of heat-shock genes is rapidly induced after temperature stress. Hsp90, by itself and/or associated with Human, Mouse chaperone complexes, is a major repressor of this gene.

Antigen:Recombinant protein of human HSF1

Form:

How to use:1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution:PBS (pH7.4) containing 1% BSA

Application:This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: