

Anti-Human/Mouse ING2 Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA06928

3th Edition

Description: This gene is a member of the inhibitor of growth (ING) family. Members of the ING family associate with and modulate the activity of histone acetyltransferase (HAT) and histone deacetylase (HDAC) complexes and function in DNA repair and apoptosis.

Antigen: Recombinant protein of human ING2

Form:

How to use: 1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution: PBS (pH7.4) containing 1% BSA

Application: This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: Widely expressed. Higher expressed in colon-cancer tumor than in normal colon tissues.