

Anti-Human/Mouse/Rat DNAJC3 Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA07027

3th Edition

Description: This gene encodes a protein with multiple tetratricopeptide repeat (TPR) motifs as well as the highly conserved J domain found in DNAJ chaperone family members. It is a member of the tetratricopeptide repeat family of proteins and acts as an inhibitor of the interferon-induced, dsRNA-activated protein kinase (PKR).

Antigen: Recombinant protein of human DNAJC3

Form:

How to use: 1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution: PBS (pH7.4) containing 1% BSA

Application: This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: Widely expressed with high level in the pancreas and testis. Also expressed in cell lines with different levels.