

**Anti-Human/Mouse SERPINF1 Polyclonal Antibody**

**Polyclonal Antibody**

**Cat.NO.: PA07046**

---

3th Edition

**Description:**The protein encoded by this gene is a member of the serpin family, although it does not display the serine protease inhibitory activity shown by many of the other serpin family members. The encoded protein is secreted and strongly inhibits angiogenesis. In addition, this protein is a neurotrophic factor involved in neuronal differentiation in retinoblastoma cells.

**Antigen:**Recombinant protein of human SERPINF1

**Form:**

**How to use:**1.0 ml distilled water will be added to the product

**Stability:** Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

**Dilution:**PBS (pH7.4) containing 1% BSA

**Application:**This antibody can be used for western blotting in concentration of 1?5?g/ml.

**Specificity:**Retinal pigment epithelial cells and blood plasma.