

**Anti-Human PRMT7 Polyclonal Antibody**

**Polyclonal Antibody**

**Cat.NO.: PA07060**

---

3th Edition

**Description:**Protein arginine methyltransferase 7 is a protein that in humans is encoded by the PRMT7 gene. Arginine methylation is an apparently irreversible protein modification catalyzed by arginine methyltransferases, such as PMT7, using S-adenosylmethionine (AdoMet) as the methyl donor. Arginine methylation is implicated in signal transduction, RNA transport, and RNA splicing. May be involved in etoposide-induced cytotoxicity, a chemotherapeutic agent frequently used for testicular cancer and small-cell lung cancer that can cause cytotoxicity in the treatment of other cancers. Down-regulation confers increased sensitivity to the Top1 inhibitor camptothecin (CPT).

**Antigen:**Recombinant protein of human PRMT7

**Form:**

**How to use:**1.0 ml distilled water will be added to the product

**Stability:** Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

**Dilution:**PBS (pH7.4) containing 1% BSA

**Application:**This antibody can be used for western blotting in concentration of 1?5?g/ml.

**Specificity:**