

Anti-Human/Mouse WNK2 Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA07132

3th Edition

Description:The protein encoded by this gene is a cytoplasmic serine-threonine kinase that belongs to the protein kinase superfamily. The protein plays an important role in the regulation of electrolyte homeostasis, cell signaling survival, and proliferation. Alternative splicing results in multiple transcript variants.

Antigen:Synthetic peptide of human WNK2

Form:

How to use:1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution:PBS (pH7.4) containing 1% BSA

Application:This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity:Expressed in various cancer cell lines (at protein level). Predominantly expressed in heart, brain, skeletal muscle and colon.