

Anti-Human VWCE Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA07148

3th Edition

Description: May be a regulatory element in the beta-catenin signaling pathway and a target for chemoprevention of hepatocellular carcinoma.

Antigen: Synthetic peptide of human VWCE

Form:

How to use: 1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution: PBS (pH7.4) containing 1% BSA

Application: This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: Expressed in liver.