

Anti-Human/Mouse/Rat VGF Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA07164

3th Edition

Description: This gene is specifically expressed in a subpopulation of neuroendocrine cells, and is upregulated by nerve growth factor. The structural organization of this gene is similar to that of the rat gene, and both the translated and the untranslated regions show a high degree of sequence similarity to the rat gene. The encoded secretory protein also shares similarities with the secretogranin/chromogranin family, however, its exact function is not known.

Antigen: Synthetic peptide of human VGF

Form:

How to use: 1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution: PBS (pH7.4) containing 1% BSA

Application: This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: Central and peripheral nervous systems, synthesized exclusively in neuronal and neuroendocrine cells.