

**Anti-Human/Mouse/Rat UCP3 Polyclonal Antibody**

**Polyclonal Antibody**

**Cat.NO.: PA07196**

---

3th Edition

**Description:** Mitochondrial uncoupling proteins (UCP) are members of the larger family of mitochondrial anion carrier proteins (MACP). UCPs separate oxidative phosphorylation from ATP synthesis with energy dissipated as heat, also referred to as the mitochondrial proton leak. UCPs facilitate the transfer of anions from the inner to the outer mitochondrial membrane and the return transfer of protons from the outer to the inner mitochondrial membrane.

**Antigen:** Synthetic peptide of human UCP3

**Form:**

**How to use:** 1.0 ml distilled water will be added to the product

**Stability:** Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

**Dilution:** PBS (pH7.4) containing 1% BSA

**Application:** This antibody can be used for western blotting in concentration of 1?5?g/ml.

**Specificity:** Only in skeletal muscle and heart. Is more expressed in glycolytic than in oxidative skeletal muscles.