

**Anti-Human HAUS7 Polyclonal Antibody**

**Polyclonal Antibody**

**Cat.NO.: PA07203**

---

3th Edition

**Description:** This gene encodes a subunit of the augmin complex, which regulates centrosome and mitotic spindle integrity, and is necessary for the completion of cytokinesis. The encoded protein was identified by interaction with ubiquitin C-terminal hydrolase 37. Alternative splicing results in multiple transcript variants.

**Antigen:** Synthetic peptide of human HAUS7

**Form:**

**How to use:** 1.0 ml distilled water will be added to the product

**Stability:** Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

**Dilution:** PBS (pH7.4) containing 1% BSA

**Application:** This antibody can be used for western blotting in concentration of 1?5?g/ml.

**Specificity:** Detected in spleen, thymus, testis, ovary, small intestine and colon, with highest levels of expression in testis and ovary.