

**Anti-Human/Mouse TSPYL1 Polyclonal Antibody**

**Polyclonal Antibody**

**Cat.NO.: PA07236**

---

3th Edition

**Description:**The protein encoded by this gene is found in the nucleolus and is similar to that of a family of genes on the Y-chromosome. This gene is intronless. Defects in this gene are a cause of sudden infant death with dysgenesis of the testes syndrome (SIDDT).

**Antigen:**Synthetic peptide of human TSPYL1

**Form:**

**How to use:**1.0 ml distilled water will be added to the product

**Stability:** Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

**Dilution:**PBS (pH7.4) containing 1% BSA

**Application:**This antibody can be used for western blotting in concentration of 1?5?g/ml.

**Specificity:**Expressed in testis, ovary, liver, spleen, brain, kidney, prostate, lung, liver, and heart.